


INTEGRATED FLUIDIC CIRCUITS

A versatile tool that concurrently automates hundreds to thousands of reactions per sample within a single run

Standard BioTools™ integrated fluidic circuits (IFCs) do the work for you. Fast-track genomic studies by eliminating upfront prep work and significantly reducing hands-on time, enabling high-complexity analysis with unprecedented proficiency. Microfluidics-based solutions using IFCs streamline workflows and substantially reduce costs for applications demanding sensitivity and broad dynamic range, including genotyping, gene expression, sample identification, copy number variation and NGS library preparation.

Proven proprietary technology: How it works

Experiments are easily transferable to IFCs and require a fraction of the sample and reagents needed for microplates. IFCs are precisely fabricated in microscopic dimensions with extremely accurate volumetric chambers for nanoscale processing of individual reactions.



Highlights

IFCs

Scalable

Multiple IFC formats make it easy to manage dynamic sample or assay counts.

Flexible

Open architecture allows the user to easily add, remove, replicate or change assays within a panel.

Cost-effective

Reduce cost and waste through lower consumption of reagents and lab plastics.

Go beyond multiplex output
with singleplex simplicity

Microfluidic arrays: Future forward

- Perform multiple PCR-based applications including targeted NGS library preparation
- Achieve lower cost and higher operating efficiency through volume miniaturization and workflow automation
- Benefit from flexible and scalable sample throughput without changing technologies
- Automation-compatible design enables integration with lab robotic systems





Plate-based PCR

Each sample (up to 96) is interrogated **ONLY** by assays within each inlet. Throughput and performance subject to multiplexing challenges.


IFC

Each sample is interrogated by **EVERY** assay in independent reaction chambers using one dye channel. Up to 96 samples per run. Up to 96 assays per sample. Up to **9,216** results from one run.



96 PLATES
VS.
1 CHIP






MICROfluidics
MEGA data
MEGA savings



IFCs are run on the Biomark™ X9 System for High-Throughput Genomics. The system contains several programs that process the IFC for thermal control and fluorescence detection, enabling different IFC formats and functions for qPCR and NGS library prep applications.




Biomark X9 System – qPCR workflow

- Flexibility and scalability to adapt to multiple projects and changing needs
- Multiple IFC formats to instantly respond to changing experimental needs

| | | Applications | Samples/Targets | Partitions | Datapoints |
|---|---|--------------|-----------------|------------|------------|
|  | Flex Six™ IFC Lowest sample throughput via six partitions that each support a matrix of 12 samples by 12 assays. Each partition can be run separately or together to support target selection and protocol development | GE, qPCR, GT | 12/12 | Six | 864 |
|  | 24.192 IFC Low sample throughput with high-content capacity, which allows labs to add assays as needed | GE, qPCR, GT | 24/192 | Zero | 4,608 |
|  | 48.48 IFC Mid-capacity to support early-phase discovery and low-volume testing sessions | GE, qPCR, GT | 48/48 | Zero | 2,304 |
|  | 96.96 IFC and 96.96 GT Preamp IFC-X Superior data-generation capacity for deep sample profiling with option for integrated preamplification | GE, qPCR, GT | 96/96 | Zero | 9,216 |
|  | 192.24 IFC High sample throughput to support peak seasonal and large studies | GE, qPCR, GT | 192/24 | Zero | 4,608 |

One system, two workflows,
endless possibilities

Biomark X9 System – NGS library prep workflow

| | Number of Samples Per Run | Number of Assay Targets Per Sample | Amplicon Length | Turnaround Time from DNA to NGS-Ready Libraries (one IFC run) | Sample Throughput in 24 Hours | Number of Samples Multiplexed Per NGS Run |
|---|-----------------------------------|--|--------------------|---|-------------------------------|--|
|  | LP 8.8.6 Up to 48 | >1,600 per partition >9,600 per IFC | 150–500 base pairs | About nine hours (includes 4.5 hours walk-away automation) | Up to 96 sample libraries | Up to 1,536 using Standard BioTools barcodes |
|  | LP 48.48 IFC Up to 48 | Up to 4,800 | 150–500 base pairs | Eight hours total (includes four hours walk-away automation) | Up to 96 sample libraries | Up to 1,536 using Standard BioTools barcodes |
|  | LP 192.24 IFC Up to 192 | Up to 2,400 | 150–500 base pairs | 9.5 hours total (includes 5.5 hours walk-away automation) | Up to 384 sample libraries | Up to 1,536 using Standard BioTools barcodes |

More than 9,000 publications referencing IFCs in genomics



5,310
publications

Gene expression

- Immune profiling
- Biomarker discovery and validation
- Pharmacodynamics
- Pathways



1,690
publications

SNP genotyping

- Pharmacogenomics
- Oncology
- Inherited disease
- Quality control



721
publications

NGS library preparation

- Pathogen detection
- Oncology
- Inherited disease

Approximate numbers based on Google Scholar searches

Save money | Save reagents | Save time

Learn more about the power of microfluidics: standardbio.com/biomark-x9

Contact us at standardbio.com/contactus

Unleashing tools to accelerate breakthroughs in human health™



FLDM-01308 Rev 03 122025
IFC Portfolio Gx brochure

For Research Use Only. Not for use in diagnostic procedures.

Patent and License Information: www.standardbio.com/legal/notices. Trademarks: www.standardbio.com/legal/trademarks. Any other trademarks are the sole property of their respective owners. ©2025 Standard BioTools Inc. All rights reserved.