

SomaLogic® Authorized Site

# Lab Installation Layout

This document describes the required equipment, dimensions, and suggested layout of a lab space processing the SomaScan® Assay Kit.




**somalogic**



# Dimensions & Space Considerations

EXAMPLE LAYOUT: 10 feet x 21 feet (3m x 6.5m)



INDICATES THIS AREA CONTAINS EQUIPMENT THAT REQUIRES POWER OUTLETS

## Tecan Fluent® 780 Automated Workstation

### Dimensions:

- Height (standard Z): 1,236 mm/48.6"
- Height (long Z): 2,301 mm/90.6"
- Width: 1,650 mm/64.96"
- Depth: 785 mm/30.9"
- Weight: 152.8 kg/ 336.8 lbs
- Shipping weight: 235.8 kg/ 519.8 lbs
- Package Dimensions: 183 cm x 113 cm x 153 cm/ 72 in. x 45 in. x 60 in.
- Electrical Requirements: 100-240V (Swapable power cord)

### Site Specifications:

#### Elevator Requirements

- Door Width Minimum: 113cm (45in.)
- Height Minimum: 153cm (60 in.)
- Length Minimum: 234cm (92in.)
- Loading Capacity Minimum: 500kg/ 1000lbs

#### Door and entryway requirements

- Width Minimum: 88cm (35in.)
- Height Minimum: 66cm (26in.)

### Recommendation:

Requires table 6 ft. or larger with additional bench space for computer. Workstation should not be exposed to direct sunlight or air flow. Workstation should be in close proximity to biohazard waste and recycling containers.



## Processing Area


### Dimensions:

6 ft. (1.8m) bench space to accommodate

- 25 °C Water Bath (recommended dimensions 13.25" x 8.25" x 9.88" / 33.6cm x 21.0cm x 25.1cm)
- Centrifuge (18.5" x 21" x 13.4" / 47.0cm x 53.3cm x 34.0cm)
- GrantBio 4-Position Shaker (15" x 15" x 5.5" / 38.1cm x 38.1cm x 14.0cm)
- Rack Thawing Station (6.5" x 4.25" x 11" / 16.5cm x 10.8 x 14.0cm)
  - Typically installed in storage shelves above processing area

### Recommendation:

Used in sample and reagent preparation, hybridization set-up, and consumable storage.



## Pipetting Station


### Dimensions:

Large, empty lab bench with space to accommodate

- Pipette Storage
- Vortexer
- Microcentrifuge

### Recommendation:

Used when pipetting sample diluent, and in auxiliary preparation and slide dropping.



## Microarray Scanning & Washing Station

### Dimensions:

6 ft. (1.8m) bench space to accommodate

- Agilent Technologies SureScan Microarray Scanner (17" x 26" x 16.5") / (43cm x 67cm x 42cm) Weight: 56.8 kg/ 125.26 lbs
- SciGene Little Dipper Slide Washer (20" x 28" x 22") Weight: 26.8 kg (59lbs)
- Additional space for cracking slides
- Requires internet access
- May require fume hood during use of ACN for processing microarray slides. Check your local Environmental Health and Safety regulations



## Hybridization Area


### Dimensions:

6 ft. (1.8m) bench space to accommodate

- Hybridization Oven roughly 40" x 20" x 20" / 102cm x 51cm x 51cm Weight: 75.0 lbs (34.0 kg)
- 37 °C Incubator (10" x 14" x 16" / 25cm x 36cm x 41cm)
- Storage for reagents, flammable liquids, and high BTU waste/Waste volume:
  - Solid Biohazardous Waste: Roughly half a 50 gal trash can per assay (2100g)
  - Liquid Biohazardous Waste: 500mL
  - Flammable Waste: 350mL

### Recommendation:

Microarray processing uses acetonitrile. Area may need vent depending on your institution's regulations.






## General Storage

### Dimensions:

Space to accommodate

- Storage for consumables including robotic tips, plates, troughs, Agilent reagents, etc.
- Desiccator for storage of unused Agilent slides (recommended)
- Ambient Kit (8" x 4" x 2-1/8") will change to 8-1/4" x 4-1/4" x 2-1/4"



## 4 °C Fridge

### Dimensions:

Space to accommodate

- 1L Assay Buffer Bottles
- 4C Kits (5-5/16" x 4-3/4" x 4-15/32" per each kit box) will change to 5-5/16" x 4-3/4" x 4-15/32"




## -20 °C Freezer

### Dimensions:

Space to accommodate

- 20 °C Kit Components (8" x 4" x 2-1/8") will change to (6" x 3-1/2" x 2-1/8")
- Additional stage for Archive Plates (recommended)




## -80 °C Freezer

### Dimensions:

Space to accommodate


- Sample and control sample storage sent with kits (4-9/32" x 3-1/32" x 2-5/8") will change to (4-9/32" x 3-1/32" x 2-5/8")



## Sink

### Recommendation:

- Sink in lab to clean reusable reagent reservoir
- Fume hood to address bio-safety concerns for sample aliquoting



## Other Lab Requirements

Access to deionized (DI) water



D0005154 Rev 3:2023-11

**SomaLogic® Authorized Site - Lab Installation Layout**

SomaLogic® SomaScan® SOMAmer® SomaSignal® and associated logos are trademarks of SomaLogic Operating Co., Inc. and any third-party trademarks used herein are the property of their respective owners.

For Research Use Only (RUO). Not intended for diagnostic or patient management purposes. SomaLogic Co, Inc. is accredited to ISO 15189:2012; ISO 27001 and is a CAP-accredited CLIA laboratory.

© 2023 SomaLogic Operating Co., Inc. | 2945 Wilderness Pl, Boulder, CO 80301 | Ph 303 625 9000 | www.somallogic.com