



SomaScan® Assay LabThread™ SLX Equipment List for Authorized Sites

SOMASCAN® ASSAY EQUIPMENT LIST





INTRODUCTION

This document lists the required and optional equipment needed to run the SomaScan Assay. The first table contains the required equipment specific to the SomaScan Assay. None of the components can be changed to different part numbers and vendors. The second table lists required laboratory equipment with vendor and part numbers listed as examples, but equivalent equipment is acceptable. The third table lists optional equipment that is used at SomaLogic and which enables higher throughput or enhanced sample tracking.

REQUIRED FLUID HANDLING SYSTEM (TABLE 1)				
VENDOR	PART NUMBER	ITEM	QUANTITY	UNDER WARRANTY
SomaLogic	420-00001	SomaLogic LabThread Automation Equipment SLX	1	Yes

REQUIRED EQUIPMENT BUNDLE (TABLE 2) This bundle is purchased from SomaLogic and must be purchased as a bundle. **PART UNDER VENDOR ITEM QUANTITY NUMBER WARRANTY**

	Tygon High Purity Tubing	1		
		Assay, Waste Trough	1	Yes
		Waste Bottle	1	
		Waste Bottle Closure	1	
		Threaded Barb Fitting for Waste Trough	3	
		Screwdriver, Ball-End Hex, 4 mm Size, 11-1/2" Overall Length	1	
	000 0000	Gravity Fed Reagent Trough	1	Yes
SomaLogic	902-00020	Thermal Adapter	7	Yes
		Copper Reagent Heat Block	1	Yes
		Deck Segment, Photocleavage Mounting	1	Yes
		Deck Segment, TMS Shaker Mounting	1	Yes
		Deck Segment, Heat Plate Mounting	3 1 1 7 1	Yes
		Slide Dropper	1	Yes
		LabThread Thermal Magnetic Shaker	3	Yes
		LabThread Thermal Magnetic Shaker Control Box	1	Yes
		Thermal Adapter Ejection Tool	1	Yes

REQUIRED EQUIPMENT BUNDLE (TABLE 2) continued

This bundle is purchased from SomaLogic and must be purchased as a bundle.

VENDOR	PART NUMBER	ITEM	QUANTITY	UNDER WARRANTY
		4 Serial Port Expansion Device	1	Yes
		USB cable for 4 serial port expansion device	1	Yes
		International Power Strip with 6 outlets	1	Yes
		AC Patch Cord - Short	3	Yes
		AC Patch Cord - Long	1	Yes
		Hook and Loop Cable Ties, 6" Overall Length, 3/8" Wide	1	
		Flame-Retardant Hook and Loop, Adhesive Backing, 1"W x 5'L	1	
		LED UV Photocleavage Lamp Mounting Bracket	1	Yes
		LED UV Lamp	2	Yes
		Ribbon Cable Assembly for LED UV Lamp	1	Yes
SomaLogic	902-00020	Disposable Biomedical Waste Container	1	
		Cardboard Recyclables Waste Container	1	
		18-8 Stainless Steel Shoulder Screw, 5/16" Shoulder Diameter	2	Yes
		18-8 Stainless Steel Flange Nut, 10-24 thread size 5/16"	2	Yes
		Cabinet Divider Wall	1	Yes
		Wireless Keyboard and Mouse Combo	1	Yes
		Monitor and Keyboard Mount	1	Yes
		Computer Monitor, 24 inch	1	Yes
		Mounting Bracket for LabThread SLX	1	Yes
		Display Port Cable 3M	1	Yes
		USB Cable A to B 3M	1	Yes
		Screw-in Hook and Loop Cable Tie Mount, 12/Pk	1	

REQUIRED EQUIPMENT BUNDLE (TABLE 2) continued

This bundle is purchased from SomaLogic and must be purchased as a bundle.

NOTE - These items are the Ambient items for the qualification kit - it wil be a separate line item in an order

VENDOR	PART NUMBER	ITEM	QUANTITY	UNDER WARRANTY
	Tube, 500 µL 2D Barcoded, V Bottom White Scr Top, Sterile, 480/cs	Tube, 500 µL 2D Barcoded, V Bottom White Screw Top, Sterile, 480/cs	1	
		Foil Plate Seal, 100/Pk	1	
		Gaskets, 8 HD Arrays/Slide	16	
		Adapter 96-Well Plate	60	
		Plate, 96 well 2 mL Deep-Well, Non-Sterile	15	
		Plate, Corning 96-Well, Black w/Clear Bottom - 1 Plate	25	
		Lids, Microplates, 100/pk	1	
		Troughs, Beckman Auto Reservoir, 96 Channel, 1 Trough	10	
		Plate, Nunc 96-Well, Round Bottom, 500 µL, 1 plate	50	
		Trough, Tecan 96-Well, Deep Well Reagent Reservoir - 1 Trough	50	
		Basin, Biologix, PVC, 100-Pack/Each	1	
SomaLogic	801-00004	Plate, Nunc 96-Well, Conical Bottom, 300µL	30	
		Filled Bottle, 2x Hybridization Buffer	4	
		SomaScan 11K Assay v5.0 - Components, Ambient 2 Slide Kit	6	
		ART 20 SoftFit-L (Rainin LTS Tips) - 1 rack	10	
		ART 200 SoftFit-L (Rainin LTS Tips) - 1 rack	30 4 6 10	
		ART 300 SoftFit-L (Rainin LTS Tips) - 1 rack	10	
		ART 1000 SoftFit-L (Rainin LTS Tips) - 1 rack	1	
		Bottle, Square PETG Media, 125mL	4	
		Bottle, Square PETG Media, 60mL	4	
		Bottle, Square PETG Media, 30mL	4	
		Tube, Screw Top w/o Cap, 2mL	2	
		Screw Cap for Microtube, w/ O-ring, Polypropylene (Yellow)	2	
		Tube, Matrix 1mL 2D Screw Top, Amber PP, V Bottom w/o Caps,	2	

REQUIRED EQUIPMENT BUNDLE (TABLE 2) continued

This bundle is purchased from SomaLogic and must be purchased as a bundle. NOTE - These items are the Ambient items for the qualification kit - it wil be a separate line item in an order

VENDOR	PART NUMBER	ITEM	QUANTITY	UNDER WARRANTY
SomaLogic		Bottle, Narrow-Mouth PETG, 5mL	4	
		Tube, 500uL 2D Barcoded, V Bottom White Screw top, Sterile	6	
	801-00004	Tube, 500uL 2D Barcoded, V Bottom Green Screw Top, Sterile	6	
		Tube, 500uL 2D Barcoded, V Bottom Red Screw Top, Sterile	10	
		Caps, Screw Tube, Clear, Sterile	2	
		Bottle, Round w/Cap, 10mL, PETG	2	

REQUIRED EQUIPMENT (TABLE 3)

This equipment is required and cannot be replaced with different vendors and/or part numbers.

VENDOR	ITEM	QUANTITY	PART NUMBER	NOTES
	SureScan Microarray Scanner	1ea	G4900DA	
Agilont	Hybridization Chamber Kit - SureHyb Enabled	12 ea	G2534A	
Agilent Technologies	Hybridization Oven	1ea	G2545A	
	Hybridization Oven Rotator Rack	1ea	G2530- 60029	
SPT Lab Tech	BioMicrolab Thaw	1ea	BioMicrolab Thaw	
Cole-Parmer	Staining Dish System; 50 Slide Capacity	1ea	EW-48585-10	
Rainin	Rainin Span 6 Pipettor (300 µL)	1ea	LA6-300XLS	
Grant Instruments	Quad Microplate Shaker-Incubator (GrantBio 4-Position Shakers)	1ea	PHMP-4	Either one GrantBio shaker or three Eppendorf Thermomixer C are needed.
Thermo Fisher Scientific	Matrix Tube Decapper	1ea	4105MAT	

	REQUIRED EQUIPMENT - SLIDE WASH OPTIONS					
VENDOR	ITEM	QUANTITY	PART NUMBER	NOTES		
SciGene (Option 1: Automatic Wash)	Little Dipper Slide Washer	1ea	1080-40-1			
	Magnetic hotplate stirrer, 5 L, with temperature probe	1ea	76549-918	Hot plate with temp control		
	5L stir plate	1ea	76549-916	Stir plate for RT washes		
	16-32 Slide Rack, Glass	1ea	53500-434	Slide Rack for 16 slide capacity (Case of 3)		
VWR (Option 2:	Dish	1ea	25461-016	Slide dish (Case of 3)		
Manual Wash)	Cover	1ea	25461-018	Dish cover for ACN (Case of 1)		
	Stir Bars	3 ea	58947-106	Stir bar x3		
	Assembly Clamp for H3770 and H3710	1ea	76640-758	Stirrer support rod clamping kit (rod, 2 clamps and hook connector		

STANDARD LABORATORY EQUIPMENT (TABLE 4)

This equipment is required but suggested parts can be replaced with equivalent equipment.

VENDOR	ITEM	QUANTITY	PART NUMBER	NOTES
Eppendorf	Centrifuge 5804, Keypad, Non-Refrigerated, with Rotor A-2-DWP, US 3	1ea	22628203	Ability to centrifuge matrix tube racks and 96 well plates. Dimension of the racks and tubes is 5.029±0.005" (L) x 3.364±0.005" (W) x 1.745±0.010" (H)
VWR	Digital Incubator (37 °C)	1ea	414005-130	Minimum interior dimensions: 10" D x 14" W x 16" H
Molecular Devices	Fluorescence and Absorbance Plate Reader	1ea	Varies	Ability to do bottom read with excitation wavelength of 535 nm and detect emission wavelength of 575 nm with a cutoff wavelength of 570 nm. Used in QC step of the eluate before hybridization. Absorbance reading capabilities required if analysis of urine samples is desired.
	Vortex Mixer	1ea		
	p2008 or 12-Channel Pipette, Manual	1ea		
	p208 or 12-Channel Pipette, Manual	1ea		
	p1000 Single Channel Pipette, Manual	1ea		
	4°C Refrigerator	1ea		
No specific vendor	-80 °C Freezer	1ea		
requirements	-20 °C Freezer	1ea		
	Flammable Cabinet	1ea		Acetonitrile storage
	Lab Timer	1ea		
	Biohazard Waste Containers	1ea		
	Microcentrifuge	1ea		Capable of spinning 1.5-2 mL screw top tubes

OPTIONAL EQUIPMENT (TABLE 5)

This optional equipment is used by SomaLogic to enable higher throughput and enhanced sample tracking.

* Cannot be replaced with different vendors and/or part numbers.

VENDOR	ITEM	QUANTITY	PART NUMBER	NOTES
Thermo Fisher	VisionMate High Speed Barcode Reader (or similar)	1ea	312800	Used in 2D barcode scanning to track sample tube samples.
Scientific	ALPS 50 V-Manual Heat Sealer	1ea	AB-1443A	Recommended for long term storage of archive plates.
Life Technologies*	Magnetic-Ring Stand (96-Well)	1ea	AM10050	Used in hybridization from archive plate when controls are hybridized for QC purposes.
Agilent Technologies	Agilent Hybridization Chamber SureHyb Enabled, Stainless	12 additional	G2534A	For use when running two assay runs per day.
U.S. Plastic Corp*	Tubing for Refillable Reagent Trough (1/4" ID x 3/8" OD x 1/16" Wall)	20 feet	57579	Used in maintenance.
Eppendorf*	ThermoMixer C Eppendorf with Eppendorf Smart Block Plates and Eppendorf Thermotop	3 ea		Either one GrantBio shaker or three Eppendorf Thermomixer C are needed. Can be replaced by Grant Instruments Quad Microplate Shaker- Incubator 4-position.
Powervar	Powervar Security II UPM	1ea	ABCE1442-11	Used as UPS battery backup, can be
Tecan	Powervar Security II UPM	1ea	30106582	purchased from either Powerwar or Tecan.
Thomas Scientific	Dry-Keeper Desiccator with Gas Ports	1ea	1226D47	Storage of unused Agilent slides once package is opened.



